

PL001-5001

NEW!
SPACIOUS
HEAVY DUTY

Patent, Innovative, Multifunction Enrich Your ♥ Life

Unique Ultra-light Weight
Power Wheelchair
Easy to Maneuver

Features

- Manufactured by ISO7176
- Ultra-light weight: N.W. 33KGs , makes portability and transportation easy and convenient.
- Supporting passenger weighting up to 150kgs
- Areospace alumium and well-constructed frames achieves the solidity and safety
- Powful and safe Lithium Battery, 12V x 2 - 20AH
- Powderful , dedicated 24V brushless motor drive
- Comfortable, Spacious, Wide Cushion

We are Offering

- Unique , excellent quality products with multiple models
- Competitive price
- Professional customer services
- Good warranty terms

Company's Certificates



Innovative & Unique Power Wheelchair



4001

5001

3004

1002

Suzhou KD Intelligent Device Co., Ltd
No. 36, Gugang Road, Chenxiang Town
Taicang City, Jiangsu Province 215400

Tel: 0512-53110088 Ext 8001
Fax: 0512-53110099
E-mail: yiwenz@yahoo.com

Ultra-light weight-
33kg Support
weighing up to **150**kg

Removable
Wide **18"** Cushion

Powerful Motor Drive
with Lithium Battery
12Vx2-20AH

Cruise mileage
21KM

Easy fold to Cars & SUV



Seat Width
450mm

12" Rear
RubberWheel

Making the Life Easy and Convenient

Innovative and Unique Power Wheelchair



- Powerful Dedicated Brushless Motor Drive,
- 24v x 250w 12" Rear Rubber Wheel

Length	Width	Height	Seat Depth	Seat Width	Seat Height	Backrest Height	N.W.
1080mm	686mm	980mm	420mm	450mm	500mm	500mm	33kg

Special designed packaging or Travel case that make it well protected during transportation

Travel Case

Special Designed Packaging

Foldable in seconds, the answer for portability

Folded Size only **440x686x785mm**

Easy to fold & extremely convenient for traveling on any type of vehicle whether it is a Car, Truck, SUV, Taxi, Bus, Train, Cruise Ship or Airplane.