

# JAY XTREME ACTIVE

Skin protection for active users



# SKIN PROTECTION AND STABILITY FOR THE ACTIVE USER

The JAY Xtreme Active cushion is a lightweight and low profile cushion for an active performance and easy transfers. Designed to provide pressure relief and stability to active users. The cushion is really easy to maintain.

The **JAY fluid pad** maximises skin integrity protection and reduces shear on the active user. The **dual cover system** promotes air flow and dissipates heat and moisture effectively. Thighs positioning can be achieved through retrofittable **leg support components**.







## LEG SUPPORT COMPONENTS

The Xtreme Active offers optional field-adjustable leg support components. The 2,5 cm high adductor and abductor components can be easily placed on the inner cover through velcros. So the thighs can be properly positioned for many clinical applications which encourages orthopedic alignment and increases the seating tolerance.

→ **Leg support kit**  
contains two adductors  
and one abductor.

Outer cover

Inner cover

Lightweight foam base



→ Microclimatic cover



→ Stretch cover



→ Incontinent cover

## DUAL-COVER TECHNOLOGY

- › The JAY Xtreme Active cushion features a dual-cover system. This innovative cover technology of an inner and outer cover provides a good microclimate and at the same time an effective incontinence management.
- › The inner cover material is water-resistant and easy to wipe off and clean, while its Aqua-guard zipper and anti-wicking thread ensure that the foam base stays dry. At the same time the microclimatic outer cover utilise 3DX™ fabric to dissipate heat and moisture to maintain the skin integrity.
- › Optionally the JAY Xtreme Active is also available with an outer stretch or incontinent cover. All cushions comes with an useful accessory bag stitched to the cover.
- › Each inner and outer cover fits for both fluid pad sizes and allow the leg supports to be fitted without creating surface tension.
- › All covers are machine washable at 60° C and quick to dry.



## THE MODULARITY OF THE XTREME ACTIVE CUSHION

The JAY Xtreme Active cushion offers two different fluid options.

**Standard fluid pad:** JAY fluid insert which offers moderate to high level of pressure redistribution for most of the users. It goes with the Standard base.

**Large fluid pad:** Designed for users with advanced muscle atrophy who do not want to go to a heavier high profile cushion but require more fluid under their bony areas. The large fluid pad has a higher fluid amount for optimal immersion and pressure redistribution of the prominent areas of the pelvic. It goes with a base with larger pelvic load area.



→ **Standard fluid pad**  
with standard base.



→ **Large fluid pad**  
in case of advanced muscle atrophy. It comes with a base with a larger pelvic load area.



## TECHNICAL DATA

<b>Cushion width</b>	34 cm – 50 cm (in 2 cm increments)
<b>Cushion depth</b>	34 cm – 50 cm (in 2 cm increments)
<b>Cushion height</b>	5 cm
<b>User weight capacity</b>	150 kg
<b>Cushion weight</b>	1,7 kg (40 x 40 with standard fluid); 2,2 kg (40 x 40 with large fluid)
<b>Inserts</b>	Standard & Large JAY Flow Fluid™
<b>Outer covers</b>	Microclimatic, stretch and incontinent
<b>Leg support kit</b>	Contains 2 adductors & 1 abductor; each 2,5 cm high (optional)

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.eu](http://www.SunriseMedical.eu) or alternatively is available on request in large text.



### Sunrise Medical GmbH

Kahlbachring 2-4  
D-69254 Malsch/HD  
Tel.: +49 (0) 72 53 / 9 80-0  
Fax: +49 (0) 72 53 / 9 80-222  
E-Mail: [international@sunrisemedical.de](mailto:international@sunrisemedical.de)  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)