

GMFCS
Levels I-III

kg
max. 60



(-)3°-17,5°



(-)5°-15°



16-40 cm



17-40 cm

I am your seating aid for everyday use.

I promote active sitting and thus support your development. The great thing is, however big you are at the moment, I will grow with you! You can simply adjust my seat and foot plate up or down or to the rear so that you can sit comfortably at the table with your big brother. I have been manufactured with loving attention to detail using wood. I suit the other furniture in your home perfectly! With my many accessory parts, I can easily be adapted even more individually to your requirements.







**Fast
height adjustment**
via the aluminium profiles
integrated in the frame.



**Fast
angular adjustment**
via pull strap





With a great
range of
accessories!



Data	Size 0 KITA	Size 0	Size 1 KITA	Size 1	Size 2
Seat width	17 – 28 cm	17 – 28 cm	20 – 35 cm	20 – 35 cm	30 – 40 cm
Seat depth	16 – 26 cm	16 – 26 cm	25 – 35 cm	25 – 35 cm	30 – 40 cm
Max. seat height	45 cm	65 cm	45 cm	65 cm	65 cm
Back height	25 cm	25 cm	29 cm	29 cm	34 cm
Back angle	(-) 3° – 17.5°				
Seat angle	(-) 5°, 0°, 7.5°, 15°				
Foot angle	(-) 5°, 0°, 7.5°, 15°				
Foot support height min.	13 cm* / 18 cm** / 21 cm***				
Foot support height min. (with fast angle adjustment)	8 cm* / 13 cm** / 16 cm***				
max. load	30 kg	30 kg	45 kg	45 kg	60 kg
Weight	9.2 kg	9.7 kg	9.6 kg	10.2 kg	12 kg
HMV No.	26.11.05.0034				

Adjustable lower leg length max. = selected/required seat height – foot support height min. * Standard, ** With 50 mm double castor, *** With 75 mm double castor



You can find more information on fabrics and embroidery
on pages 174 – 175.

You can find more information on accessories and harnesses
in the Price List.

