

Push handle
for easy transfer.

Back unit
optionally removable, can be
adjusted in angle and height.

Foot support
for good positioning of the feet.

Castors
available in different sizes,
facilitate transfer.

Headrest
can be adjusted in height,
depth and angle.

Lateral pads
allow a stable and
upright sitting position.



Therapy table
offers a lot of space
for activities.

Seat surface
is individually adaptable.

Foot pedal
for height adjustment.

Base frame
compact and sturdily built.



**In just one single step, my
back unit can simply be
removed. This gives you
more freedom to move.**



Frame colours

● Anthracite metallic

● Telegrey 4



Data	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
Back height with short back	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
Back height with long back	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
Seat height	25 - 62 cm	25 - 62 cm	25 - 62 cm	25 - 62 cm
Seat height with electric height adjustment	32 - 62 cm	32 - 62 cm	32 - 62 cm	32 - 62 cm
Seat height with hydraulic height adjustment	27 - 63 cm	27 - 63 cm	27 - 63 cm	27 - 63 cm
Lower leg length	17 - 24 cm	22 - 32 cm	27 - 38 cm	30 - 48 cm
Total width	54 cm	54 cm	64 cm	64 cm
Length of base frame	61 cm	61 cm	75 cm	75 cm
Total length with combination sub-frame	76 cm	76 cm	86 cm	86 cm
Seat angle	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°
Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
max. load	30 kg	30 kg	45 kg	50 kg
Weight	19 kg	20 kg	22 kg	24 kg
HMV No.	26.11.05.0022			



You can find more information on fabrics & embroidery
on pages 174 - 175.

You can find more information on accessories & harnesses
in the Price List.

