

GMFCS
Levels III-V

kg
max. 50



(-)5°–25°

(-)5°–37°

15–42 cm

46–58 cm

16–38 cm



I am the therapy chair that grows with you.

I really love accompanying you as a small child or adolescent. That's why I am available in four different sizes. The great thing about me is that my seat can tilt forwards or to the rear – it's really easy to do this using a lever. And I can do this so far that you can actually lie down and relax. I can be raised and lowered so that you can always be in the centre of whatever is going on. I can also help you to train your back musculature: to do this, remove my backrest – it's a quick and easy process.







**In just one single step:
if you like, my back unit
can simply be removed.
This gives you more free-
dom to move.**



Frame colours

- Anthracite metallic
- Telegrey 4

**Removable
back unit**



Data	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 – 25 cm	21 – 30 cm	26 – 36 cm	30 – 42 cm
Seat width	16 – 26 cm	20 – 30 cm	20 – 30 cm	28 – 38 cm
Back height with short back	35 – 43 cm	35 – 43 cm	42 – 52 cm	42 – 52 cm
Back height with long back	42 – 52 cm	42 – 52 cm	53 – 63 cm	53 – 63 cm
Seat height	46 – 58 cm	46 – 58 cm	46 – 58 cm	46 – 58 cm
Seat height with low sub-frame	35 – 40 cm	35 – 40 cm	-	-
Lower leg length	17 – 24 cm	22 – 32 cm	27 – 38 cm	30 – 48 cm
Total width	45 cm	45 cm	53 cm	53 cm
Length of base frame	66 cm	66 cm	77 cm	77 cm
Seat angle	(-)5° – 37°	(-)5° – 37°	(-)5° – 37°	(-)5° – 37°
Back angle	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°
max. load	30 kg	30 kg	45 kg	50 kg
Weight	16 kg	17 kg	18 kg	19 kg
HMV No.	26.11.05.0021			



You can find more information on fabrics and embroidery
on pages 174 – 175.

You can find more information on accessories and harnesses
in the Price List.

