

activity chair.

Simply moving.

GMFCS
Levels II-IV



I am your dynamic therapy chair.

I really hate sitting still! I'd much rather be changed into lots of different positions. I can tilt forwards or backwards, up or down – as dynamically as you require. Active sitting has a major advantage: it helps you to relieve physical and psychological tension. My wide range of accessories such as headrests, pelotte pads or side guides provide the necessary stability in any situation.





activity chair.

Push handle
for convenient transfer.

Headrest
can be adjusted in depth, height
and angle.

Thorax pelotte pads
facilitate a stable and up-
right sitting position.

Back angle
can be adjusted using a
gas pressure spring.

Armrests
can be adjusted in height and
angle.

Thigh guide
can be adjusted in angle.

Abduction wedge
guides the legs.

Seat pad
provides more comfort.

Sub-frame
is particularly sturdy
and adjustable in height.

Foot support
can be folded up, and is
therefore ideal for indepen-
dent transfer.

Dual castors with locking devices
facilitate transfer.

Optional dynamic function
gently absorbs your movements
and can be easily switched on
and off using the rotary handle
without requiring tools.



Frame colour ● Effect black



Data	Size 0	Size 1	Size 2
Seat depth	20 - 30 cm	28 - 41 cm	38 - 51 cm
Seat width	30 cm	36 cm	43 cm
Seat width (with hip pelotte pads)	18 - 23 cm	22 - 29 cm	28 - 36 cm
Chest width	14 - 29 cm	17 - 30 cm	24 - 36 cm
Back height	32 - 39 cm	39 - 48 cm	48 - 61 cm
Headrest height - seat	37 - 53 cm	44 - 62 cm	50 - 75 cm
Armrest height - seat	13 - 19 cm	18 - 27 cm	18 - 27 cm
Lower leg length	23 - 30 cm	30 - 41 cm	38 - 51 cm
Seat height over the floor With short / with long feet	24 - 32 cm / 47 - 55 cm	32 - 42 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
Seat height over the floor With short / with long feet and dual castors	34 - 42 cm / 47 - 55 cm	41 - 51 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
Frame width With short / with long feet	53 cm / 58 cm	59 cm / 62 cm	66 cm / 66 cm
Frame width With short / with long feet and dual castors	57 cm / 60 cm	62 cm / 65 cm	71 cm / 71 cm
Seat angle	(-)15° - 15°	(-)15° - 15°	(-)15° - 15°
Back angle	(-)5° - 20°	(-)5° - 20°	(-)5° - 20°
max. load	34 kg	68 kg	113 kg
HMV No.	26.11.05.0xxx		



You can find more information on fabrics and embroidery
on pages 174 - 175.

You can find more information on accessories and har-
nesses
in the Price List.

