



SALSA R²
COMPACT AND
POWERFUL





> Easy transfer

The smart flip-up armrest design enables easy transfer. It can be infinitely adjusted in height, depth and angle.

// FULL RANGE OF POWERED SEATING AND SPECIAL CONTROL OPTIONS

Change your seating position. You can do this with three different types of tilt modules – the powered or mechanical single pivot tilt option for a basic tilt, the powered tilt module with integrated centre of gravity adjustment and the lift & tilt combi module.

The backrest can be optional reclined with an actuator or a mechanical gas strut.

Finally the full range of special controls or input devices can be added to the chair to find the best operating performance for the user.

FULLY ADJUSTABLE SEATING SYSTEM

The Salsa R² is equipped with a fully adjustable seating system which can be quickly adjusted throughout the full range of width and depth, back angle and height without additional parts. The seat is designed for perfect adaptation to the user and enables prompt servicing and refurbishment. Alternative backs and cushions can be easily fitted.



Powered tilt module
(0° to 30°) with centre of
gravity adjustment to suit
correct users balance.





New base integrated LED lights and indicators (optional)

4 point tie-down system for transportation in a vehicle. Optional with Dahl car docking system available.

Adjustable drive wheel suspension



- **Stylish and compact base**
in 3 colours available - blue, red & white

// **SALSA R² IS APPROVED FOR TRANSPORT IN A VEHICLE**

The Salsa R² is transit tested according to ISO 7176-19 and therefore approved for transportation in a vehicle.

The tests are done for the 4 point tie-down system and the optional Dahl docking system for an even easier fixation.

ENHANCES YOUR MOBILITY. INSIDE AND OUTSIDE.



> Active Swing Away hangers

Swing away hangers for easy transfer. Also available as elevating legrest (manual or powered)

The Salsa R2 shows a great indoor and outdoor performance. Indoors the compact base and the low seat to floor height starting from 42cm ensures easy under-table access. Outdoors the adjustable suspension and the powerful motors (2-pole or 4-pole) always offer a comfortable and active driving characteristic up to a speed of 10 kph. With the simple battery access from the rear the batteries can be easily changed while the user is in the chair.

Simple battery and controller access while user is in the chair

Batteries available with 60 or 70 Ah

12,5" and 14" drive wheels

9" and 10" castor wheels



SALSA R²

Compact and powerful



Configure the wheelchair to your individual needs - The Salsa R² offers a wide range of options.
For further information on product options please check the Salsa R² order form.



**// POWERED TILT MODULE WITH
CENTRE OF GRAVITY ADJUST-
MENT**

offers a tilt from 0° up to 30° with an additional weight shift mechanism to suit correct users balance.



**// POWERED SINGLE PIVOT TILT
OPTION**

for basic tilt from 3° to 20° to quickly change between a driving and a relax seat position



// POWERED RECLINE BACK

This recline backrest offers a range from -5° to 30°. Easily retrofittable if developing needs requires higher support.



**// POWERED ARTICULATING
LEGREST**

for advanced leg support. Also available as powered elevating legrest



// TRUNK SUPPORT LATERALS

Multi-adjustable laterals offer advanced positioning of the trunk.



**// FULL RANGE OF CONTROLS
DEVICES AVAILABLE**

- VR2 standard control (up to two actuators)
- Up to four actuator control through a separate switch box
- Dual control
- R-Net optional controller
- Full range of special controls

TECHNICAL DATA

SALSA R² // power wheelchair

Max. user weight:	140 kg	Battery size:	60 and 70Ah
Caster wheel options:	9" and 10"	Max. dynamic stability:	8° (14%)
Drive wheel options:	12" and 14"	Overall width:	60 cm (12" wheels); 62 cm (14" wheels)
Seat width:	41 – 51 cm	Turning radius:	max. 120 cm
Seat depth:	41 – 51 cm	Max. kerb climbing:	10 cm (with kerb climber)
Max. speed:	6 kph and 10 kph	Colour choices:	3 colours - red, blue & white
Max. range (ISO 7176-4):	32 km (40 km with 70 Ah)	Transport (ISO 7176-19):	4 Point tie-down & Dahl docking system

For further information on product options please check the order form on our website.
You can view this document also in a larger font on a PDF-file at www.SunriseMedical.eu



Sunrise Medical GmbH

Kahlbachring 2-4
D-69254 Malsch/Heidelberg
T +49 (0) 7253 980-0
F +49 (0) 7253 980-222
www.sunrisemedical.eu

