

Design



COMFORT

- New Ergonomic Design
- Suspension Structure



SAFE

- Strong AL6061-T6 Welding Frame
- Hydro-Foaming Main-Tube
- Intelligent Electromagnetic Brake

SGS FDA CE CERAH



CONTROL

- Brushless Control System
- Easy Operation
- Soft Start & Stop
- Foldable & Adjustable Joystick



BATTERY

- More than 24.3Km Cruise mileage
- 24V12AH Ternary lithium Battery
- More than 800 Charging Cycle times
- 10+ Protection Circuit board
- UN38.3 certificate
- MSDS report
- Air travel friendly



POWER DRIVE

- 2-Power Drive
- Foldable Anti-tipping wheels
- Load up to 120kgs



FOLD

- Fold and unfolding in seconds
- Easy to operate
- Small compact size
- Easy storage and transportation

Use for



Park



Travel



Shopping



Working



Indoor

Features



19.8kgs Without Battery



Pull To Fold



Pull To Unfold



Pull For External Charge



Switch To Manual Mode



Lift To Front On



Press To Side On



Foldable & Adjustable Joystick



Trolley Mode



7001

Weight /Without Battery
22.3kg/19.8kg
L x W x H (Unfolded)
870*620*910mm (34.2"*24.4"*35.8")
L x W x H (Folded)
800*620*320mm (31.4"*24.4"*12.6")
User Weight
120kg
Rear Wheel
8"
Front Wheel
6"
Battery
12AH
Package Size
675*475*890mm (26.6"*18.7"*35")
G.W.
22.3g (19.8kg)



7002

Weight /Without Battery
23kg/20.5kg
L x W x H (Unfolded)
900*620*930mm (35.4"*24.4"*34")
L x W x H (Folded)
800*620*320mm (31.4"*24.4"*12.6")
User Weight
120kg
Rear Wheel
10"
Front Wheel
7"
Battery
12AH
Package Size
675*475*890mm (26.6"*18.7"*35")
G.W.
23kg (20.5kg)



7003

Weight /Without Battery
25.4kg/22.9kg
L x W x H (Unfolded)
930*640*950mm (36.6"*25.1"*37.4")
L x W x H (Folded)
800*640*320mm (31.4"*25.4"*12.6")
User Weight
120kg
Rear Wheel
12"
Front Wheel
7"
Battery
12AH
Package Size
720*495*890mm (28.3"*19.5"*35")
G.W.
25.4kg(22.9kg)



Suzhou KD Intelligent Device Co., Ltd
No. 36, Gugang Road, Chenxiang Town Taicang City,
Jiangsu Province 215400

Tel: 0512-53110088 Ext 8001
Fax: 0512-53110099
E-mail: yiwenz@yahoo.com

The Pursue of Mobility...